

Protocol: Yeast DAPI Stain

DAPI is a fluorescent nucleic acid stain. In this protocol ethanol is used as a fixative.

What is needed:

A yeast culture to stain
Phosphate buffered saline (PBS)
Stock DAPI solution (1mg/ml)
Agarose pads for microscopy and coverslips

1. Centrifuge 1ml of yeast for 1 minute at 6,000 x g.
2. Decant supernatant. Wash in 1ml of PBS. Centrifuge for 1 minute at 6,000 x g.
3. Decant supernatant and resuspend in 1ml PBS.
4. Add 333 μ l of the washed yeast culture to a new 1.5ml microcentrifuge tube.
5. Add 666 μ l of 100% ethanol to the yeast culture.
6. Incubate the yeast and ethanol mixture at room temperature for 30 minutes. (You may notice clumping of the cells; this is OK as ethanol causes proteins to congeal or stick together.)
7. While the cells are being fixed, make the DAPI solution by adding 1.25 μ l of a 1mg/ml DAPI stock solution to 1ml of phosphate buffered saline (PBS). Keep this solution in the dark until it is needed in **Step 10** below.
8. Centrifuge the cells for 1 minute at 2,500 rpm.
9. Decant the supernatant. Wash the cells in 1ml of PBS without DAPI. Centrifuge again for 1 minute at 2,500 rpm.
10. Decant the supernatant. Add 200 μ l of the DAPI stain and resuspend the cells.
11. Add 2 μ l (up to 5 μ l depending on cell density) of cells to a microscope slide with agarose pad. Add a drop of VectaShield (if desired). Place coverslip and seal with nail polish.